

<https://www.iqniter.com/release-notes/>

IQNITER Studio v3.9.1 (Cardio Training)

New Features

- Added 'Switch Monitors' button to Settings panel enabling Group Screen and Instructor Screen to be switched.
- Added in the Session tab option for the instructor to specify additional Music Search Folder.

Optimizations and bug-fixes

- Fixed that RPM was not always shown correctly on the Group Screen.
- Fixed that exerciser tiles did not always show correct values and color at end of session.
- Optimized switching layout to/from Session Profile layout (toggle button F3) so the session timers do not freeze for a few seconds.
- Added Incline to Cue Card.

IQNITER DesiQner v3.9.1

Optimizations and bug-fixes

- Session music can now be saved on USB devices again.
- Now the Search field is in focus when the Session Manager is opened.